

PREMIUM RESTORATIVES

QUADRANT UNIVERSAL LC

Quadrant Universal Living Colors (QULC) is a genuine all-round composite; a universal hybrid suited to multiple types of layering techniques.

Thanks to its balanced composition, QULC can be applied and processed easily, and lends itself perfectly to major restorations. What makes this product unique is the fluorescent pigment, which simulates tooth structure, ensuring that QULC complements the natural and vibrant appearance of the dentition. Thanks to the superior quality of the pigments and monomers used, QULC guarantees long-lasting colour stability.

Available in the following shades:
A1, A2, A3, A3.5, A4, OA2, B2, B3, C2, C3, I.



PREMIUM RESTORATIVES

YOUR IMPRESSION IS OUR SPECIALTY



ZYF200.01

PREMIUM RESTORATIVES

QUADRANT

The Cavex Quadrant system delivers ideal filling materials for all situations. Experience the proven and reliable quality.



INNOVATIVE RESTORATIVE SOLUTIONS FOR PROFESSIONAL USE

Discover the Quadrant system; your gateway to a smooth and efficient treatment yielding amazing results! That is why Quadrant is the pre-eminent choice for dentists dedicated to delivering extraordinary restoration work.

Infinite possibilities

The Quadrant system encompasses a wide range of colours and viscosities, giving dentists the freedom to carefully select the appropriate filling for each individual patient. Whether for aesthetic restorations in the visible region, or for functional repairs further back in the mouth, the Quadrant system is sure to provide the right solution. In addition, all Quadrant composites and bonding systems are chemically compatible, offering excellent options for mutual combination.

Consistency and durability

Our Quadrant stands out, not only for the extensive choice it offers, but also by virtue of the consistency and durability of the materials used. Dentists can rely on the proven performance of the products, while patients benefit from high quality and long-lasting restorations in terms of both function and aesthetics. Did you know that Quadrant composites feature a unique fluorescent pigment? This is an optical equivalent of true dental tooth structure, giving restorations a natural and vibrant appearance.

Ease of use as a priority

The Quadrant system provides a solution for all dentists, and has been developed with ease of use in mind. Quadrant composites are structured as so-called hybrid composites. Average particle size is 0.7 µm, and the product is enriched with a mix of fine, pre-polymerised particles. This formula enables smooth application of the composite, as well as very easy processing, offering perfect wall adaptation. Quadrant composites can be easily modelled, do not stick to the instruments and can be polished to an impressively high gloss. Dentists can therefore work efficiently, minimising chair time for their patients.

Quality at the right price

Quadrant is available to all dentists striving for quality and perfection in their restorative work. The tried, tested and trusted quality has proved itself time and time again. At Cavex we believe that high-quality dental care should be available for everyone. With our Quadrant system, we furnish dentists with a product of exceptional quality that goes hand in hand with affordability. Our materials enable dentists to deliver compromise-free care without exceeding their budgets.

QUADRANT SYSTEM



Source: Dr. H. de Kloet

BIARIUM GLASS

The filler particles in Quadrant are prism-shaped, so that incoming light is absorbed and dispersed in all directions. As a consequence, the restoration will blend fully with the element, rendering colour differences less critical. The result: natural and invisible restorations. Thanks to the optimal light dispersion, only one colour composite is required in most cases (mono layer technique).

QUADRANT ANTERIOR SHINE

Quadrant Anterior Shine is the right choice for those professionals specialising in restorations of the anterior part of the dentition.

This light-curing, micro-hybride composite gives your patient a healthy and natural appearance. Quadrant Anterior Shine provides gloss, polishes easily, and thanks to its natural fluorescence is barely distinguishable from natural teeth.

A special tool has been developed for the Anterior Shine composite: the 'Quadrant Quast'. This is a brush developed for the accurate application and distribution of Quadrant Anterior Shine.

Available in the following shades: A1, A2, A3, A3.5, A4, B, B2, B3, C2, C4, DS, D3 OA2, OA35, I.

QUADRANT FLOW

Quadrant Flow is not just any composite; it is distinctive thanks to its thixotropic character.

Flowing effortlessly but maintaining its shape when required, Quadrant Flow delivers the perfect balance between ease and precision. The composite filling is user-friendly, and thanks to the small syringes and their tips, the composite can be applied directly.

QUADRANT BONDING: OPTIMAL RESTORATIONS

Correct and reliable adhesion is indispensable when creating the very best restorations in dentition.

Quadrant Uni-SE-Bond is a 7th-generation, self-etching bonding system that can be combined with the Quadrant composites. By virtue of the combination of a proven 4-META monomer and mild etching, this is a strong and durable product. It is also very user-friendly, creating superb results using a simple procedure.

Quadrant Uni-1-Bond, combined with Quadrant Total Etch, is the best choice for those professionals who insist on maximum enamel adhesion and perfect visual control of the etching process. The use of 35% phosphoric acid gel for the total-etch technique offers the ideal preparation, while the 4-META in the combined primer and sealer ensures optimum bonding strength.

Cavex also offers an easy solution for the application of primers, bonding products and desensitizers: Q-brushes. These disposable brushes are ideal for one-time use.



Source: Dr Hein de Kloet

If you would like to know more about any of these products, please visit our website.

