


PREMIUM ALGINATE SYSTEM

CAVEX SCANSPRAY

Cavex ScanSpray is developed to optimize the scannability of impressions. Because it is a powder based spray it can be used on any type of impression material. Due to the addition of titanium dioxide even the smallest details are clearly scannable.



Apply a thin layer of Cavex ScanSpray to the impression for a 100% accurate digital model.



The Cavex ScanSpray sprays over 420 seconds.



ZAT01E.01

YOUR IMPRESSION IS OUR SPECIALTY

CAVEX

Fustweg 5 | 2031 CJ Haarlem | The Netherlands | cavex.nl