


Cavex ScanSpray


Cavex ScanSpray is developed to optimize the scannability of impressions. Because it is a powder based spray it can be used on any type of impression material. Due to the addition of titanium dioxide even the smallest details are clearly scannable.



Apply a thin layer of Cavex ScanSpray to the impression for a 100% accurate digital model.



The Cavex ScanSpray sprays over 420 seconds.



- For a super detailed digital model of the impression.
- Suitable for all types of impression material.
- Powder-based instead of oil-based, so it does not affect the impression.